



Końcówki Micro Tips Premium 2,5 µl, bezbarwne



Pasują do pipet:

Biozym (PreCision 2,5 µl, PreCision 10 µl),
Biohit (Proline 2,5 µl, Proline 10 µl, mLine 3 µl, mLine 10 µl),
Brand (Transferpette 1 µl, Transferpette 10 µl),
Eppendorf (2,5 µl (antracytowa), 10 µl),
Finn (Digital Micro 2 µl, Digital Micro 10 µl, Focus Micro 10 µl),
Gilson (Pipetman P-2, Pipetman P-10)
Nichiryo (2 µl, 10 µl),
Rainin (Pipet-Classic PR2, Pipet-Lite SL-2).
Wolne od DNaz i RNaz.

Miniatura	Nr-art.	Nazwa	Pojemność	Kolor	Stopień czystości	Pakowane
	A-720025	Końcówki Micro Tips Premium	2,5 µl	bezbarwne	niesterylne	w worku (1000 szt.)
	A-720028	Końcówki Micro Tips Premium	2,5 µl	bezbarwne	niesterylne	w statywie (10 statywów po 96 szt.)
	A-720029	Końcówki Micro Tips Premium	2,5 µl	bezbarwne	sterylne	w statywie (10 statywów po 96 szt.)

Opis

Pasują do pipet:

Biozym (PreCision 2,5 µl, PreCision 10 µl),
Biohit (Proline 2,5 µl, Proline 10 µl, mLine 3 µl, mLine 10 µl),
Brand (Transferpette 1 µl, Transferpette 10 µl),
Eppendorf (2,5 µl (antracytowa), 10 µl),
Finn (Digital Micro 2 µl, Digital Micro 10 µl, Focus Micro 10 µl),
Gilson (Pipetman P-2, Pipetman P-10)
Nichiryo (2 µl, 10 µl),
Rainin (Pipet-Classic PR2, Pipet-Lite SL-2).
Wolne od DNaz i RNaz.