

	Analogowe	Cyfrowe
Orbitalne	3005 (nośność 8 kg) → K-1590  3015 (nośność 15 kg) → K-1596  3019 (nośność 30 kg) → K-1600 	3017 (nośność 15 kg) → K-1598  3020 (nośność 30 kg) → K-1601 

	Analogowe	Cyfrowe
Posuwisto- zwrotne	3006 (nośność 8 kg) → K-1591 	3018 (nośność 15 kg) → K-1599 
	3016 (nośność 15 kg) → K-1597 	
Obrotowo- kołyszące	3011 (nośność 15 kg) → K-1592 	3012 (nośność 15 kg) → K-1593 
	3013 (nośność 15 kg) → K-1594 	3014 (nośność 15 kg) → K-1595 
Kołyszące		

	Analogowe	Cyfrowe
Kołowo-wibracyjne	3023 (nośność 1,2 kg – na płytki mikrotitracyjne) → K-1602 	
Obrotowe	3040 (nośność 20 kg – na butelki) → K-1604  3025 (max 24 probówki Ø 12-17 mm, dł. 75-180 mm) → K-1603 	